


Product Datasheet

SULT1C2 antibody (orb1336466)

Description	Mouse monoclonal antibody to SULT1C2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human SULT1C2 (NP_001047) produced in HEK293T cell.
Target	SULT1C2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	34.7 kDa
Uniprot ID	O00338
NCBI	NM_001056
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...