


Product Datasheet

HDHD1A (PUDP) antibody (orb1336558)

Description	Mouse monoclonal antibody to HDHD1A (PUDP)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human HDHD1(NP_036212) produced in HEK293T cell.
Target	PUDP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	23.5 kDa
Uniprot ID	Q08623
NCBI	NM_012080
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...