


Product Datasheet

MAGEB18 antibody (orb1336613)

Description	Mouse monoclonal antibody to MAGEB18
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human MAGEB18 (NP_775970) produced in HEK293T cell.
Target	MAGEB18
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	38.3 kDa
Uniprot ID	Q96M61
NCBI	NM_173699
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...