


Product Datasheet

VASP antibody (orb1336631)

Description	Mouse monoclonal antibody to VASP
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human VASP (NP_003361) produced in HEK293T cell.
Target	VASP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	39.6 kDa
Uniprot ID	P50552
NCBI	NM_003370
Dilution Range	WB 1:500~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...