


Product Datasheet

EIF4EL3 (EIF4E2) antibody (orb1336667)

Description	Mouse monoclonal antibody to EIF4EL3 (EIF4E2)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human EIF4E2 (NP_004837) produced in HEK293T cell.
Target	EIF4E2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	28.2 kDa
Uniprot ID	O60573
NCBI	NM_004846
Dilution Range	WB 1:1000~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...