


Product Datasheet

RTF2 antibody (orb1336674)

Description	Mouse monoclonal antibody to RTF2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human C20orf43 (NP_057491) produced in HEK293T cell.
Target	RTF2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	33.7 kDa
Uniprot ID	Q9BY42
NCBI	NM_016407
Dilution Range	WB 1:1000~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...