

Product Datasheet

PKC epsilon (PRKCE) antibody (orb1336679)

Description Mouse monoclonal antibody to PKC epsilon (PRKCE)

Species/Host Mouse

Reactivity Human, Monkey, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full length human recombinant protein of human PRKCE (NP_005391) produced in HEK293T cell.

Target PRKCE

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2b

Clonality Monoclonal


MW 83.5 kDa

Uniprot ID [Q02156](#)

NCBI [NM_005400](#)

Dilution Range WB 1:1000~2000, IHC 1:150, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...