


Product Datasheet

VCAM1 antibody (orb1336730)

Description	Mouse monoclonal antibody to VCAM1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human VCAM1(NP_001069) produced in HEK293T cell.
Target	VCAM1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	78.7 kDa
Uniprot ID	P19320
NCBI	NM_001078
Dilution Range	WB 1:500~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...