


Product Datasheet

RDH14 antibody (orb1336773)

Description	Mouse monoclonal antibody to RDH14
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RDH14 (NP_065956) produced in HEK293T cell.
Target	RDH14
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	36.7 kDa
Uniprot ID	Q9HBH5
NCBI	NM_020905
Dilution Range	WB 1:200~500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...