


Product Datasheet

Butyrylcholinesterase (BCHE) antibody (orb1336779)

Description	Mouse monoclonal antibody to Butyrylcholinesterase (BCHE)
Species/Host	Mouse
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human BCHE (NP_000046) produced in HEK293T cell.
Target	BCHE
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	68.2 kDa
Uniprot ID	P06276
NCBI	NM_000055
Dilution Range	WB 1:500~2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...