


Product Datasheet

CTDSP1 antibody (orb1336795)

Description	Mouse monoclonal antibody to CTDSP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CTDSP1 (NP_067021) produced in E.coli.
Target	CTDSP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	29 kDa
Uniprot ID	Q9GZU7
NCBI	NM_021198
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...