


Product Datasheet

Gephyrin (GPHN) antibody (orb1336829)

Description	Mouse monoclonal antibody to Gephyrin (GPHN)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human GPHN (NP_065857) produced in HEK293T cell.
Target	GPHN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	83.3 kDa
Uniprot ID	Q9NQX3
NCBI	NM_020806
Dilution Range	WB 1:1000~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...