


Product Datasheet

PACSIN3 antibody (orb1336836)

Description	Mouse monoclonal antibody to PACSIN3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human PACSIN3 (NP_057307) produced in HEK293T cell.
Target	PACSIN3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	48.3 kDa
Uniprot ID	Q9UKS6
NCBI	NM_016223
Dilution Range	WB 1:200~500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...