


Product Datasheet

SQSTM1 antibody (orb1336838)

Description	Mouse monoclonal antibody to SQSTM1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human SQSTM1 (NP_003891) produced in HEK293T cell.
Target	SQSTM1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	47.5 kDa
Uniprot ID	Q13501
NCBI	NM_003900
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...