


Product Datasheet

NMNAT1 antibody (orb1336855)

Description	Mouse monoclonal antibody to NMNAT1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human NMNAT1 (NP_073624) produced in HEK293T cell.
Target	NMNAT1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	31.8 kDa
Uniprot ID	Q9HAN9
NCBI	NM_022787
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...