

Product Datasheet

EXDL1 (EXD1) antibody (orb1336912)

Description	Mouse monoclonal antibody to EXDL1 (EXD1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human EXD1 (NP_689809) produced in HEK293T cell.
Target	EXD1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	58.2 kDa
Uniprot ID	Q8NHP7
NCBI	NM_152596
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.

