


Product Datasheet

HDAC1 antibody (orb1336926)

Description	Mouse monoclonal antibody to HDAC1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human HDAC1 (NP_004955) produced in HEK293T cell.
Target	HDAC1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	54.9 kDa
Uniprot ID	Q13547
NCBI	NM_004964
Dilution Range	WB 1:200~500, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...