


Product Datasheet

CNDP1 antibody (orb1336970)

Description	Mouse monoclonal antibody to CNDP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CNDP1 (NP_116038) produced in HEK293T cell.
Target	CNDP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	56.5 kDa
Uniprot ID	Q96KN2
NCBI	NM_032649
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...