


Product Datasheet

SLFN1 antibody (orb1336982)

Description	Mouse monoclonal antibody to SLFN1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SLFN1 (NP_659427) produced in HEK293T cell.
Target	SLFN1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	45.4 kDa
Uniprot ID	Q499Z3
NCBI	NM_144990
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...