


Product Datasheet

PUS7 antibody (orb1337011)

Description	Mouse monoclonal antibody to PUS7
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PUS7 (NP_061915) produced in HEK293T cell.
Target	PUS7
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	74.9 kDa
Uniprot ID	Q96PZ0
NCBI	NM_019042
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...