


Product Datasheet

PBLD antibody (orb1337166)

Description	Mouse monoclonal antibody to PBLD
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PBLD (NP_071412) produced in HEK293T cell.
Target	PBLD
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	31.6 kDa
Uniprot ID	P30039
NCBI	NM_022129
Dilution Range	WB 1:2000, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...