


Product Datasheet

SULT1A1 antibody (orb1337207)

Description	Mouse monoclonal antibody to SULT1A1
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human SULT1A1 (NP_001046) produced in HEK293T cell.
Target	SULT1A1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	34.0 kDa
Uniprot ID	P50225
NCBI	NM_001055
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...