


Product Datasheet

DNAJA2 antibody (orb1337279)

Description	Mouse monoclonal antibody to DNAJA2
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human DNAJA2 (NP_005871) produced in HEK293T cell.
Target	DNAJA2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	45.6 kDa
Uniprot ID	O60884
NCBI	NM_005880
Dilution Range	WB 1:500~2000, IHC 1:50~150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...