


Product Datasheet

PIK3AP1 antibody (orb1337294)

Description	Mouse monoclonal antibody to PIK3AP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human PIK3AP1 (NP_689552) produced in HEK293T cell.
Target	PIK3AP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	90.2 kDa
Uniprot ID	Q6ZUJ8
NCBI	NM_152309
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...