


Product Datasheet

NSE2 (NSMCE2) antibody (orb1337369)

Description	Mouse monoclonal antibody to NSE2 (NSMCE2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human NSMCE2 (NP_775956) produced in E.coli.
Target	NSMCE2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	27.8 kDa
Uniprot ID	Q96MF7
NCBI	NM_173685
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...