


Product Datasheet

FSH beta (FSHB) antibody (orb1337378)

Description	Mouse monoclonal antibody to FSH beta (FSHB)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human FSHB (NP_001018090) produced in HEK293T cell.
Target	FSHB
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	12.4 kDa
Uniprot ID	P01225
NCBI	NM_001018080
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...