


Product Datasheet

CBWD1 antibody (orb1337393)

Description	Mouse monoclonal antibody to CBWD1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CBWD1 (NP_060961) produced in HEK293T cell.
Target	CBWD1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	Q9BRT8
NCBI	NM_018491
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...