


Product Datasheet

mu Crystallin (CRYM) antibody (orb1337504)

Description	Mouse monoclonal antibody to mu Crystallin (CRYM)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CRYM (NP_001879) produced in HEK293T cell.
Target	CRYM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	33.6 kDa
Uniprot ID	Q14894
NCBI	NM_001888
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...