


Product Datasheet

Ephrin A2 (EFNA2) antibody (orb1337527)

Description	Mouse monoclonal antibody to Ephrin A2 (EFNA2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human EFNA2 (NP_001396) produced in HEK293T cell.
Target	EFNA2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	21.3 kDa
Uniprot ID	O43921
NCBI	NM_001405
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...