

Product Datasheet

IL36 alpha (IL36A) antibody (orb1337555)

Description Mouse monoclonal antibody to IL36 alpha (IL36A)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human IL1F6 (NP_055255) produced in HEK293T cell.

Target IL36A

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal

MW 17.5 kDa

Uniprot ID [Q9UHA7](#)

NCBI [NM_014440](#)

Dilution Range WB 1:2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.

