

Product Datasheet

CFAP298 antibody (orb1337632)

Description	Mouse monoclonal antibody to CFAP298
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human C21orf59 (NP_067077) produced in HEK293T cell.
Target	CFAP298
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	33.2 kDa
Uniprot ID	P57076
NCBI	NM_021254
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...