


Product Datasheet

SH3PX1 (SNX9) antibody (orb1337653)

Description	Mouse monoclonal antibody to SH3PX1 (SNX9)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human SNX9(NP_057308) produced in HEK293T cell.
Target	SNX9
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	66.4 kDa
Uniprot ID	Q9Y5X1
NCBI	NM_016224
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...