


Product Datasheet

HIBCH antibody (orb1337657)

Description	Mouse monoclonal antibody to HIBCH
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human HIBCH (NP_055177) produced in HEK293T cell.
Target	HIBCH
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	39.4 kDa
Uniprot ID	Q6NVY1
NCBI	NM_014362
Dilution Range	WB 1:1000~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...