


Product Datasheet

HIBCH antibody (orb1337660)

Description	Mouse monoclonal antibody to HIBCH
Species/Host	Mouse
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human HIBCH (NP_055177) produced in HEK293T cell.
Target	HIBCH
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	39.4 kDa
Uniprot ID	Q6NVY1
NCBI	NM_014362
Dilution Range	WB 1:1000~2000, IHC 1:50, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...