


Product Datasheet

IGF2BP2 antibody (orb1337702)

Description	Mouse monoclonal antibody to IGF2BP2
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human IGF2BP2 (NP_006539) produced in HEK293T cell.
Target	IGF2BP2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	65.9 kDa
Uniprot ID	Q9Y6M1
NCBI	NM_006548
Dilution Range	WB 1:500~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...