


Product Datasheet

PTP epsilon (PTPRE) antibody (orb1337740)

Description	Mouse monoclonal antibody to PTP epsilon (PTPRE)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Immunogen	Full length human recombinant protein of human PTPRE (NP_006495) produced in HEK293T cell.
Target	PTPRE
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	80.5 kDa
Uniprot ID	P23469
NCBI	NM_006504
Dilution Range	IHC 1:150, IF 1:100, FLOW 1:50~100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...