


## Product Datasheet

### DAP Kinase 2 (DAPK2) antibody (orb1337765)

|                            |                                                                                                                                                                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to DAP Kinase 2 (DAPK2)                                                                                                                                                                                                    |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | FC, IF, IP, WB                                                                                                                                                                                                                                       |
| <b>Immunogen</b>           | Full length human recombinant protein of human DAPK2 (NP_055141) produced in HEK293T cell.                                                                                                                                                           |
| <b>Target</b>              | DAPK2                                                                                                                                                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                |
| <b>Application notes</b>   | For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. |
| <b>Isotype</b>             | IgG2b                                                                                                                                                                                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                                                                           |
| <b>MW</b>                  | 42.7 kDa                                                                                                                                                                                                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q9UIK4</a>                                                                                                                                                                                                                               |
| <b>NCBI</b>                | <a href="#">NM_014326</a>                                                                                                                                                                                                                            |
| <b>Dilution Range</b>      | WB 1:2000, IF 1:100, FLOW 1:100, IP 2ug/500ul                                                                                                                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                      |



immunofluorescent staining of COS7 cells...