


Product Datasheet

Glucose 6 phosphate isomerase (GPI) antibody (orb1337795)

Description	Mouse monoclonal antibody to Glucose 6 phosphate isomerase (GPI)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human GPI(NP_000166) produced in HEK293T cell.
Target	GPI
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	63.0 kDa
Uniprot ID	P06744
NCBI	NM_000175
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...