

Product Datasheet

Cytochrome P450 1A2 (CYP1A2) antibody (orb1337809)

Description Mouse monoclonal antibody to Cytochrome P450 1A2 (CYP1A2)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, WB

Immunogen Full length human recombinant protein of human CYP1A2 (NP_000752) produced in HEK293T cell.

Target CYP1A2

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 58.2 kDa

Uniprot ID [P05177](#)

NCBI [NM_000761](#)

Dilution Range WB 1:500, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...