


Product Datasheet

NME4 antibody (orb1337825)

Description	Mouse monoclonal antibody to NME4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human NME4 (NP_005000) produced in HEK293T cell.
Target	NME4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.6 kDa
Uniprot ID	O00746
NCBI	NM_005009
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...