


Product Datasheet

TUBA8 antibody (orb1337866)

Description	Mouse monoclonal antibody to TUBA8
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human TUBA8(NP_061816) produced in HEK293T cell.
Target	TUBA8
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q9NY65
NCBI	NM_018943
Dilution Range	WB 1:2000, IHC: 1:50-1:150, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...