


Product Datasheet

PIM2 antibody (orb1337892)

Description	Mouse monoclonal antibody to PIM2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human PIM2 (NP_006866) produced in HEK293T cell.
Target	PIM2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG3
Clonality	Monoclonal
MW	34 kDa
Uniprot ID	Q9P1W9
NCBI	NM_006875
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100, IP: 4ug/mL
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...