


Product Datasheet

TUBA8 antibody (orb1337924)

Description	Mouse monoclonal antibody to TUBA8
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	Full length human recombinant protein of human TUBA8(NP_061816) produced in HEK293T cell.
Target	TUBA8
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q9NY65
NCBI	NM_018943
Dilution Range	WB 1:2000; IF 1:100; IP: 4ug/mL
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...