

Product Datasheet

PTP epsilon (PTPRE) antibody (orb1337935)

Description Mouse monoclonal antibody to PTP epsilon (PTPRE)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human PTPRE (NP_006495) produced in HEK293T cell.

Target PTPRE

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal

MW 80.6 kDa

Uniprot ID [P23469](#)

NCBI [NM_006504](#)

Dilution Range WB 1:500, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.

