


Product Datasheet

MTRF1L antibody (orb1337957)

Description	Mouse monoclonal antibody to MTRF1L
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human MTRF1L (NP_061914) produced in HEK293T cell.
Target	MTRF1L
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	43.6 kDa
Uniprot ID	Q9UGC7
NCBI	NM_019041
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...