


Product Datasheet

LIPG antibody (orb1337963)

Description	Mouse monoclonal antibody to LIPG
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human LIPG (NP_006024) produced in HEK293T cell.
Target	LIPG
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	54.6 kDa
Uniprot ID	Q9Y5X9
NCBI	NM_006033
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...