


Product Datasheet

KIAA1609 (TLDC1) antibody (orb1337998)

Description	Mouse monoclonal antibody to KIAA1609 (TLDC1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human KIAA1609 (NP_065998) produced in HEK293T cell.
Target	MEAK7
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	51.0 kDa
Uniprot ID	Q6P9B6
NCBI	NM_020947
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...