


Product Datasheet

Annexin A11 (ANXA11) antibody (orb1338009)

Description	Mouse monoclonal antibody to Annexin A11 (ANXA11)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human ANXA11 (NP_001148) produced in HEK293T cell.
Target	ANXA11
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	54.4 kDa
Uniprot ID	P50995
NCBI	NM_001157
Dilution Range	IHC 1:100~200, WB: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...