


Product Datasheet

Pleckstrin (PLEK) antibody (orb1338057)

Description	Mouse monoclonal antibody to Pleckstrin (PLEK)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human PLEK (NP_002655) produced in HEK293T cell.
Target	PLEK
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	40.1 kDa
Uniprot ID	P08567
NCBI	NM_002664
Dilution Range	WB 1:2000, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...