


Product Datasheet

RALBP1 antibody (orb1338059)

Description	Mouse monoclonal antibody to RALBP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Full length human recombinant protein of human RALBP1 (NP_006779) produced in HEK293T cell.
Target	RALBP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	76.0 kDa
Uniprot ID	Q15311
NCBI	NM_006788
Dilution Range	WB 1:2000, Flow 1:100, IP: 4ug/mL
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...